

TECHNICAL SPECIFICATIONS

COMFORT VALUE Range
Tilt-in-Space (CVTIS)

ISSUE DATE: 11-10-19
VERSION: TS/CVTIS-01
STATUS: APPROVED
AUTHORISED BY: Technical Director



K.YOUNG



CONTENTS

1. Application
2. Features Overview
3. Product Dimensions
4. Warranty Information
5. Product Parts List
6. Other supporting documentation

This Technical Specification has been issued by Care & Independence Ltd. It is intended in relation to the product stated only. Care & Independence reserve the right to amend the specifications at any given time. All data is correct at time of issue date. In the event of any doubt, please contact us.

Care & Independence Ltd Registered office no 4214056. T. 01869 327999 W. careandindependence.com

©Care&IndependenceLtd 2019 All Rights Reserved.

1. APPLICATION













The COMFORT VALUE TIS chair is designed to be used as attendant operated.

2. FEATURES OVERVIEW

FEATURES (CVTIS)	
Height Adjustable Frame	3" adjustment (76mm)
Wheel configuration	4x125mm (5") dia. braked castors
Seat type	Removable / Reversible wipe clean, vinyl upholstered comfort gap seat
Back Rest	Tension adjustable wipe clean vinyl
Foot Rests	Height adjustable/swing away/ folding wheelchair style
Head Rest	Fully removable wipe clean vinyl
Arm Rests	Rise & Fall movement, self-supporting fully padded with hand grips
Tilt Mechanism	Gas strut with solid rod (no cables). 45° travel.
Frame Material	25mm dia. Stainless Steel (SS304) tubular frame, white powder coated
Fixtures & Fittings	All either fully galvanised / stainless steel materials
Colour	Black Seat, Grey backrest, Black footrest plates
Manufacture	Fabricated & Assembled in the UK by Care & Independence
Warranty	12 months
Wash Care	All parts can be wiped clean using mild soap solution / cleaning agent. Do not use bleach. Do not immerse. Wipe dry after washing. Suitable for decontamination and up to 80°C wash temps.

ACCESSORIES (GENERAL)	
Lap strap	Up to 44" waist as standard (longer by quote). Black safety-belt material with side squeeze plastic clip.
Commode Pan/Lid/ Rack	3.5L capacity, oval shaped, unisex design vessel complete with anti-spill lid and fittings.
Foot Strap	Footplate mounted webbing straps with hook & loop adjustment.
Thoracic Support	Wheelchair style harness w/ waist belt, shoulder straps, groin strap to assist in maintaining an upright seated position.
Lateral Head Support	Removable wedge to fit headrest
Extended Warranty	Available on request by quote

3. PRODUCT DIMENSIONS

PRODUCT DIMENSIONS (CVTIS)											
Seat width	Overall width	Seat height	Overall height	Seat depth	Overall depth	Overall depth when seated	Armrest height	Handle clearance	Headrest height only	SWL	Weight
											
20" (508 mm)	620 mm	510 mm lowest	1275mm lowest	500 mm	1095 mm	1350 mm	180 mm	980 mm	250 mm	180kg	17kg

4. WARRANTY INFORMATION

This product is offered with a standard 12 months warranty from the date of purchase, subject to the following general terms and conditions:

- i. Should problems be experienced with any workmanship or materials within the first 12 months please contact your point of purchase or care and Independence.
- ii. Please quote the serial number / Batch of the product (where applicable), which will be found on the front cover on this user guide or on the product itself.
- iii. The product will be repaired free of charge provided that it has been used in accordance with the guidelines and recommendations set out in the user guide, and maintained by Care & Independence or a competent person
- iv. The warranty does not cover replacements; adjustments, repairs arising from a consequence of normal wear & tear, wilful damage, accidental damage, misuse, neglect or any cause, which is beyond the manufacturers control.
- v. The warranty does not cover modifications or repairs made to the product by unauthorised persons or the use of unauthorised parts or any failure or malfunction due to incorrect and/or poor workmanship on the part of any person, other than Care and independence or a competent person
- vi. This warranty does not affect your statutory consumer rights.

5. PRODUCT PARTS LIST (CVTIS)

ACCESSORIES (CVTIS)	
TITLE	PART NO
Commode Pan & Rack	PCOMPAN
Replacement Rack	PCOMRACK
Black Webbing Lap strap	LAPSTRAP/BK
Thoracic Support	ACLATSUP48B
Foot Strap	FOOTSTRAP

SPARE PARTS (CVTIS)	
TITLE	PART NO
Headrest Cover	CVTIS/002
Backrest	CVTIS/001
Calf strap	CVTIS/003
Castor - single unit	CV/001
E Clip – single	PASCLIP
Gas Strut	CAMLOCK/V1
Gas Strut Rubber Boot	RCAMBOOT
Rubber arm Pad	RBKARMPAD
Gap Front Seat - 19"	CVGF/19
Hand Grip	R1BHG
Footrest -Pair	CV/002
Commode Pan Rack – pair	PCOMRACK

6. OTHER SUPPORTING DOCUMENTATION

- User Guide - available as a download from www.careandindependence.com
- Pricing - available from the relevant Care & Independence Area Manager
- Test Results - available upon written request from Care & Independence Ltd